

SMS[®]

From Himalayas to the world!!

MANUAL KERATOMETER



Internal Scale Reading

Features easy to align mires with internal readings in one field of view simplifying and speeding measurements



LED Lamp Design

Designed for minimum maintenance, enhanced brightness, Long Life and less heating



Joystick

Unlike traditional keratometer, SMK-8586 facilitates precise and fast measurement with joystick system



Contact Lens Measurement

Able to measure radius of curvature of a contact lens with contact lens holder

**2 YEARS
WARRANTY**

SAFE MEDICAL SUPPLY

Balaju, Kathmandu, Nepal
Email: smsnepalpl@gmail.com



नेपालमा बनेको

SMS[®]

From Himalayas to the world!!

MANUAL KERATOMETER

SMK-8586



SPECIFICATIONS MODEL: SMK-8586

Scale	Corneal Radius of Curvature	5.5 -12mm
	Minimum Reading	0.01mm
	Corneal refractive Power	28D-60D
	Minimum Reading	0.125D
	Axis of corneal astigmatism	0 to 1800
	Minimum Reading	10
	Dioptre adjusting scale	-5D ~ +5D
Voltage/Frequency	220 V/50HZ	
Dimensions	330(L) X 280(W) X430(H) mm	
Weight(without motorized table)	15 KG	

SAFE MEDICAL SUPPLY

Balaju, Kathmandu, Nepal
Email: smsnepalpl@gmail.com



नेपालमा बनेको