

## SMS 2-STEP SLIT-LAMP MICROSCOPE

### SL-200 H

Binocular 2-Step Magnification

Slit Lamp Microscope

- Superior optical quality, provides stereoscopic and comfortable view, the image is clear.
- With broad beam, wide field and high depth of field.
- Lighting with Halogen light source;
- Interface for applanation tonometer reserved.



Specifications	
<b>Microscope type</b>	Cross-Angle Stereovision
<b>Magnification changer</b>	2-Step magnifier
<b>Magnification</b>	10x, 16x
<b>Eye-piece Magnification</b>	10x
<b>Visual Angle of Eye-piece</b>	13°
<b>PD Range</b>	55mm-82mm
<b>Diopter</b>	±5D
<b>Slit Projection Magnification</b>	2/3x
<b>Visual Field Diameter</b>	10x (18mm), 16x (14.5mm)
<b>Illumination Type</b>	Tower Type (Hagg-type)
<b>Light Source</b>	Halogen bulb
<b>Illumination Tilting</b>	5°, 10°, 15°, 20°
<b>Slit Width</b>	Continuous from 0mm to 9mm
<b>Slit Height</b>	Continuous from 1mm to 8mm
<b>Slit Angle</b>	0 -180° continuously adjustable
<b>Slit Aperture Size</b>	(Ø0.2, Ø1, Ø2, Ø3, Ø5, Ø8, Ø9) mm
<b>Filter</b>	Heat-Absorption, Grey, Red-free(green), Blue
<b>Voltage / Frequency</b>	115/230V ~ 50/60 Hz