

SMS[®]

DIGITAL VISION CHART MODEL: DVC-110




SAFE MEDICAL SUPPLY

BALAJU, KATHMANDU, NEPAL

Contact: smsnepalpl@gmail.com

DIGITAL VISION CHART MODEL: DVC-110



Product Features

- 21.5 inches screen, using narrow edge drawing process for frame , thin and smart
- PLS LCD panel, 30% light transmittivity, abandon frosted screen visual experience.
- High definition: visual symbol matching to complete pixel unit, completely eliminate half pixel raw edge.
- Fast response time of 0.6s, bringing smooth control experience.

Technical Parameters

Panel type	21.5-inch HD PLS panel
Resolution	1920×1080
Maximum brightness	320cd/m ²
Working distance	Adjustable 1.0-6.5m (step length 0.1m)
Installation mode	Wall mounted and Desktop Type Optional
Basic Sight mark	C Mark /E Mark /Letters/Number /ETDRS visual symbol
Functional mark	Red and Green/ Astigmatism/ Binocular balance/View Fusion, etc. Total 15 kinds
Dimension & Weight	Main frame of Wall hanging type: 519mm(W)×327.7mm(H)x70mm(L) / 4kg
	Main frame of Desktop Type: 519mm(W)x391.7mm(H)x167.5mm(L) / 3.55kg
	Remote control Panel: 55mm(W)x186mm(H) x17.2mm(L) / 80g
Power supply parameter	Input 100V-240V,50Hz-60Hz; Output DC12V,45W

SAFE MEDICAL SUPPLY

BALAJU, KATHMANDU, NEPAL

Contact: smsnepalpl@gmail.com