




SMS[®]


OPERATING MICROSCOPE

OP VISION 20



SALIENT FEATURES:

- Crisp Apochromatic Optics.
- X-Y Coupling standard.
- Easy control of X-Y position, illumination and focus level by robust foot pedal.
- Standard 0-180 Degree inclinable binocular for comfort of surgeons.
- Ergonomic handgrip for precise Handling.
- Stereo Coobervation(Optional).
- CCD HD Camera(Optional).



OPERATING MICROSCOPE

OP VISION 20



SPECIFICATIONS:

- 5 Step Manual Magnification changer.
- Motorized focus range 40mm
- 12.5X wide-field eyepieces.
- PD adjustment range: 55mm~75mm.
- 0-180 degree inclinable binocular tube with PD adjustment knob.
- Objective lens, $f=200\text{mm}$ with fine focusing.
- Magnification range: 0.4X, 0.6X, 1.0X, 1.6X, 2.5X.
- X-Y coupling range: 40mmX40mm, one button reversal to main position.

SMS[®]

SAFE MEDICAL SUPPLY

BALAJU, KATHMANDU, NEPAL

Contact: smsnepalpl@gmail.com